

<b>Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data</b>	<b>Document No.:</b>	<b>MCCRD-SOP0064</b>
	<b>Version:</b>	<b>1.0</b>
	<b>Effective Date:</b>	<b>02/03/2023</b>
	<b>Page No.:</b>	<b>Page 1 of 2</b>

### SIGNATURE APPROVALS:

Laboratory Director: Must be signed by a minimum of one of those listed at the right:	Biswajit Das, Ph.D. Chris Karlovich, Ph.D. Mickey Williams, Ph.D.	Signature / Date:
Author/Owner:	Lily Chen, Ph.D.	Signature / Date:

### REVISION HISTORY:

Document No.	Version	Description of Revision	Effective Date
MCCRD-SOP0064	1.0	Original Release	02/03/2023

### 1.0 PURPOSE/SCOPE

This Standing Operating Procedure (SOP) describes procedures for detection of Human Leukocyte Antigen (HLA) typing using whole exome sequencing (WES) data for reporting in the NCI Patient-Derived Models database as performed by the Molecular Characterization Laboratory (MoCha) at the Frederick National Laboratory for Cancer Research.

**This SOP is for research purposes only and no clinical samples will be processed using this SOP.**

### 2.0 REFERENCE DOCUMENTATION

Number	Title/Link
[1]	Szolek, A, Schubert, B, Mohr, C, Sturm, M, Feldhahn, M, and Kohlbacher, O (2014). OptiType: precision HLA typing from next-generation sequencing data Bioinformatics, 30(23):3310-6.
[2]	<a href="https://github.com/FRED-2/OptiType">https://github.com/FRED-2/OptiType</a>

### 3.0 RELATED SOPS

Available on the PDMR website: <https://pdmr.cancer.gov/sops>

Document Number	Title
MCCRD-SOP0011	Whole Exome Sequencing Data Analysis Pipeline and Specifications

**Approved For Use**

<b>Detection of Human Leukocyte Antigen (HLA) Typing from Whole Exome Sequencing Data</b>	<b>Document No.:</b>	<b>MCCRD-SOP0064</b>
	<b>Version:</b>	<b>1.0</b>
	<b>Effective Date:</b>	<b>02/03/2023</b>
	<b>Page No.:</b>	<b>Page 2 of 2</b>

#### 4.0 DESCRIPTION OF HLA Typing DETECTION

- 4.1 The mouse reads are first filtered out from whole exome sequence (WES) data following the WES data analysis pipeline in the SOP MCCRD\_SOP0011.
- 4.2 HLA reads are extracted from paired-end human only sequencing data.
- 4.3 HLA type is estimated using OptiType<sup>[1,2]</sup>.

#### 5.0 CODE DESCRIPTION

- 5.1 HLA reads are extracted from paired-end sequencing data
 

```
razers3 -i 95 -m 1 -dr 0 -tc 20 -o ${sample}_1.bam
$OPTITYPE_HOME/data/hla_reference_dna.fasta ${sample}_hg19_R1.fastq.gz
razers3 -i 95 -m 1 -dr 0 -tc 20 -o ${sample}_2.bam
$OPTITYPE_HOME/data/hla_reference_dna.fasta ${sample}_hg19_R2.fastq.gz
samtools bam2fq ${sample}_1.bam >${sample}_1.fastq
samtools bam2fq ${sample}_2.bam >${sample}_2.fastq
```

- 5.2 OptiType is used to detect HLA type for PDX sample

```
OptiTypePipeline.py -i {sample}_1.fastq {sample}_2.fastq --dna -v -o
${sample} -c $OPTITYPE_HOME/config.ini
```

- 5.3 HLA type is reported at model level with #concordance/#total PDX specimens in the model. The data for PDC and PDOrg is reported separately.

#### 6.0 Example HLA Profile Reporting in the NCI PDMR database:

- 6.1 All PDX tumors sequenced (n=5) are in concordance

HLA Profile(s)

A1	A2	B1	B2	C1	C2	PDX Concordance / PDM Type
A*02:01	A*25:01	B*08:01	B*14:02	C*07:01	C*08:02	5/5

- 6.2 In a small fraction of models, intra-model HLA variation in sequenced PDX tumors is observed.

HLA Profile(s)

A1	A2	B1	B2	C1	C2	PDX Concordance / PDM Type
A*01:01	A*32:01	B*08:01	B*14:01	C*05:01	C*07:01	1/6
A*01:01	A*32:01	B*08:01	B*14:01	C*07:01	C*08:02	3/6
A*01:01	A*32:01	B*08:01	B*14:01	C*07:01	C*08:04	2/6