

|  |                                   |             |
|--|-----------------------------------|-------------|
| SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures |                                   |             |
| Laboratory:  | Patient-Derived Models Repository |             |
| Effective Date:  | 9/6/2018                          | Page 1 of 4 |

## **Table of Contents**

|   |   |
|---|---|
| CHANGE HISTORY.....                             | 1 |
| RELATED SOPS.....                               | 1 |
| 1.0 PURPOSE/SCOPE .....                         | 2 |
| 2.0 SAFETY .....                                | 2 |
| 3.0 CLEAN-UP.....                               | 2 |
| 4.0 EQUIPMENT .....                             | 2 |
| 5.0 PREPARE 25% MATRIGEL® WORKING SOLUTION..... | 3 |
| 6.0 PREPARATION OF MATRIGEL®-COATED FLASKS..... | 3 |

## **CHANGE HISTORY**

| Revision   | Description   |
|------------|---|
|            | Internal SOP used by PDMR In Vitro Laboratory   |
| 10/15/2017 | Standardize SOP for posting to PDMR SharePoint site for use by designated NCI intramural laboratories |
| 5/14/2018  | Updated reference SOPs and Purpose/Scope section  |
| 9/6/2018   | Clarify steps in Matrigel coating and length of time for storage before use.                          |

## **RELATED SOPS**

|  |
|--|
| SOP30103: Initial Culture, Sub-culture, and Cryopreservation of Adherent Patient-Derived Tumor Cultures (PDCs) |
| SOP30105: Initial Culture and Sub-culture of Patient-Derived Cancer-Associated Fibroblasts (CAFs)              |

|  |                                   |             |
|--|-----------------------------------|-------------|
| SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures |                                   |             |
| Laboratory:  | Patient-Derived Models Repository |             |
| Effective Date:  | 9/6/2018                          | Page 2 of 4 |

## 1.0 PURPOSE/SCOPE

This Standing Operating Procedure (SOP) describes preparation of Matrigel-coated plates for successful thawing and early culture of adherent Patient-Derived Tumor Cultures (PDCs) and Cancer-Associated Fibroblasts (CAFs) under BSL-2 safety criteria. This SOP is used/performed by the Biological Testing Branch (BTB) at NCI-Frederick, Frederick National Laboratory for Cancer Research.

## 2.0 SAFETY

BTB treats all patient-derived in vitro cell cultures under Biosafety Level 2 (BSL2) conditions even when PCR-based screening has not detected the presence of a known set of human pathogens. All work is conducted in a biological safety cabinet (BSC) using personal protective equipment and avoiding the use of sharps where possible. All materials potentially exposed to the cell cultures are disinfected by exposure to a 10% bleach solution for a minimum of 10 minutes, double bagging for autoclaving or incineration. Consult with your facility safety professionals regarding the safe handling of BSL2 studies.

## 3.0 CLEAN-UP

- 3.1 All materials coming into contact with patient tissue as well as the mice carrying patient tumor samples are treated as a potential health threat (BSL-2 precautions) since the human tissues could retain human pathogenic agents even if they do not replicate in mouse cells (e.g., HIV, HPV, etc).
- 3.2 Flush/soak any items (e.g., tubes, syringes, petri dishes, lab mats, etc) that were in contact with human tissue with disinfectant (e.g., 10% bleach, commercial hydrogen peroxide disinfectant, 2% Virkon®) for a minimum of 10 minutes before disposal in biohazard waste or sharps containers (follow institutional guidelines and manufacturer's recommendations).
- 3.3 For items that can't be rinsed (e.g., micropipettors), wipe down thoroughly with bleach-soaked gauze or other appropriate disinfectants.

## 4.0 EQUIPMENT

- 4.1 Equipment
  - 4.1.1 50-mL, 25-mL, 10-mL, 5-mL pipettes, sterile
  - 4.1.2 15 and 50-mL polypropylene tubes, sterile
  - 4.1.3 Tissue Culture flasks, sterile
  - 4.1.4 Pipetman and sterile tips
  - 4.1.5 Waste container Bleach (Clorox, 5.25% Hypochlorite) diluted 1:10, 2% Virkon®, or similar disinfectant
  - 4.1.6 Refrigerator (4°C) and freezer (-20°C)
  - 4.1.7 37°C Incubator (5% CO<sub>2</sub>, humidified)
  - 4.1.8 Biological Safety Cabinet (BSC) meeting biosafety level 2 (BSL2) standards
  - 4.1.9 Personal Protective Equipment (PPE) at a minimum laboratory coat, with fitted sleeves, latex or nitrile gloves and safety glasses

|  |                                   |             |
|--|-----------------------------------|-------------|
| SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures |                                   |             |
| Laboratory:  | Patient-Derived Models Repository |             |
| Effective Date:  | 9/6/2018                          | Page 3 of 4 |

## 5.0 PREPARE 25% MATRIGEL® WORKING SOLUTION

### 5.1 Reagents

**5.1.1** 1X F-12 Nutrient Mix, without supplementation (Invitrogen, Cat#: 11765-054)

**5.1.2** Matrigel®, High concentration (BD Biosciences, Cat#: 354248)

**IMPORTANT:** All Matrigel® purchases should be submitted specifying PCR-tested LDEV-Negative Matrigel®. If not, there is a possibility of LDEV contamination which can result in LDEV+ tumors as a consequence

### 5.2 Prepare 25% Matrigel® Working Solution

**5.2.1** Chill pipettes and conical tubes in a -70°C freezer overnight and then place on wet ice prior to use.

**5.2.2** Thaw Matrigel® overnight by placing a vial of Matrigel, buried” in ice, in the refrigerator.

**5.2.3** Using cold pipettes and tubes, make a 25% Matrigel® Working Solution with 1X F-12 Nutrient Mix, without supplementation.

**5.2.4** Aliquot 25% Matrigel® either as 5-mL aliquots into chilled sterile 15-mL conical tubes or 1-mL aliquots into 2-mL sterile screw-capped tubes and place into a -20°C freezer.

## 6.0 PREPARATION OF MATRIGEL®-COATED FLASKS

### 6.1 Reagents

**6.1.1** Complete Media – PDC-specific per Certificate of Analysis (COA)

**6.1.2** 1X F-12 Nutrient Mix, without supplementation (Invitrogen, Cat#: 11765-054)

**6.1.3** Matrigel®, High concentration (BD Biosciences, Cat#: 354248)

**IMPORTANT:** All Matrigel® purchases should be submitted specifying PCR-tested LDEV-Negative Matrigel®. If not, there is a possibility of LDEV contamination which can result in LDEV+ tumors as a consequence

### 6.2 Coat Tissue Culture Plates/Flasks with 2.5% Matrigel® Solution

**6.2.1** Matrigel® coated plates/flasks should be prepared at least 1 hour before use.

- Flasks can be prepared up to 4-5 days before use and stored at 4°C following the polymerization step. Bring to ambient or 37°C before adding cells.

**6.2.2** The recommended volume of 2.5% Matrigel needed per well/flask are:

| Plate/Flask site | Volume 2.5% Matrigel <sup>®</sup> /well or flask |
|------------------|--|
| 96-well plate    | 75 µL/well                                       |
| 24-well plate    | 0.3 mL/well                                      |
| 6-well plate     | 2 mL/well  |
| T25 flask        | 3.5 mL   |
| T75 flask        | 6.0 mL   |
| T162 flask       | 8.0 mL   |

**6.2.3** The day before coating, remove the appropriate number of 25% Matrigel<sup>®</sup> Working Solution aliquots from the freezer and thaw overnight buried in ice, in the refrigerator.

**6.2.4** Prepare sufficient 2.5% solution of Matrigel<sup>®</sup> in the appropriate Complete Media (per COA) for the cells being cultured to coat vessel.

**6.2.5** Using sterile procedures, coat the growth surface of the plate/flask using the recommended volume based on the size of the well/flask (Section 6.2.2).

Be sure to rock the plate back and forth to completely coat the surface of the plate/flask.

**6.2.6** Incubate plates/flasks for a minimum of 30 minutes at ambient temperature for polymerization. This process can be enhanced by incubating in a 37°C incubator.

**6.2.7** Immediately before use, remove excess media from the flask/plate and discard taking care to not dislodge the Matrigel<sup>®</sup> coating.

- Flasks can be prepared up to 4-5 days before use and stored at 4°C following the polymerization step. Bring to ambient or 37°C before adding cells.