

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures		
Laboratory:	Patient-Derived Models Repository	
Effective Date:	1/16/2019	Page 1 of 4

Table of Contents

CHANGE HISTORY.....	1
RELATED SOPS.....	1
1.0 PURPOSE/SCOPE	2
2.0 SAFETY	2
3.0 CLEAN-UP.....	2
4.0 EQUIPMENT	2
5.0 PREPARE 25% MATRIGEL [®] WORKING SOLUTION.....	3
6.0 PREPARATION OF MATRIGEL [®] -COATED FLASKS.....	4

CHANGE HISTORY

Revision	Description
	Internal SOP used by PDMR In Vitro Laboratory
10/15/2017	Standardize SOP for posting to PDMR SharePoint site for use by designated NCI intramural laboratories
5/14/2018	Updated reference SOPs and Purpose/Scope section
9/6/2018	Clarify steps in Matrigel coating and length of time for storage before use.
1/16/2019	Added the need for Pen/Strep in the coating solution. Streamlined protocol for readability.

RELATED SOPS

SOP30103: Initial Culture, Sub-culture, and Cryopreservation of Adherent Patient-Derived Tumor Cultures (PDCs)

SOP30105: Initial Culture and Sub-culture of Patient-Derived Cancer-Associated Fibroblasts (CAFs)

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures		
Laboratory:	Patient-Derived Models Repository	
Effective Date:	1/16/2019	Page 2 of 4

1.0 PURPOSE/SCOPE

This Standing Operating Procedure (SOP) describes preparation of Matrigel-coated plates for successful thawing and early culture of adherent Patient-Derived Tumor Cultures (PDCs) and Cancer-Associated Fibroblasts (CAFs) under BSL-2 safety criteria. This SOP is used/performed by the Biological Testing Branch (BTB) at NCI-Frederick, Frederick National Laboratory for Cancer Research.

2.0 SAFETY

BTB treats all patient-derived in vitro cell cultures under Biosafety Level 2 (BSL2) conditions even when PCR-based screening has not detected the presence of a known set of human pathogens. All work is conducted in a biological safety cabinet (BSC) using personal protective equipment and avoiding the use of sharps where possible. All materials potentially exposed to the cell cultures are disinfected by exposure to a 10% bleach solution for a minimum of 10 minutes, double bagging for autoclaving or incineration. Consult with your facility safety professionals regarding the safe handling of BSL2 studies.

3.0 CLEAN-UP

- 3.1 All materials coming into contact with patient tissue as well as the mice carrying patient tumor samples are treated as a potential health threat (BSL-2 precautions) since the human tissues could retain human pathogenic agents even if they do not replicate in mouse cells (e.g., EBV, HPV, etc).
- 3.2 Flush/soak any items (e.g., tubes, syringes, petri dishes, lab mats, etc) that were in contact with human tissue with disinfectant (e.g., 10% bleach, commercial hydrogen peroxide disinfectant, 2% Virkon®) for a minimum of 10 minutes before disposal in biohazard waste or sharps containers (follow institutional guidelines and manufacturer's recommendations).
- 3.3 For items that can't be rinsed (e.g., micropipettors), wipe down thoroughly with bleach-soaked gauze or other appropriate disinfectants.

4.0 EQUIPMENT

- 4.1 Equipment
 - 4.1.1 50-mL, 25-mL, 10-mL, 5-mL pipettes, sterile
 - 4.1.2 15 and 50-mL polypropylene tubes, sterile
 - 4.1.3 Tissue Culture flasks, sterile, vented
 - 4.1.4 Pipetman and sterile tips
 - 4.1.5 Waste container Bleach (Clorox, 5.25% Hypochlorite) diluted 1:10, 2% Virkon®, or similar disinfectant
 - 4.1.6 Refrigerator (4°C) and freezer (-20°C)
 - 4.1.7 37°C Incubator (5% CO₂, humidified)
 - 4.1.8 Biological Safety Cabinet (BSC) meeting biosafety level 2 (BSL2) standards
 - 4.1.9 Personal Protective Equipment (PPE) at a minimum laboratory coat, with fitted sleeves, latex or nitrile gloves and safety glasses

SOP30102: Preparation of Matrigel-Coated Flasks for Adherent Patient-Derived In Vitro Cultures		
Laboratory:	Patient-Derived Models Repository	
Effective Date:	1/16/2019	Page 3 of 4

5.0 PREPARE 25% MATRIGEL® WORKING SOLUTION

5.1 Reagents

5.1.1 1X F-12 Nutrient Mix, without supplementation (Invitrogen, Cat#: 11765-054) with 100 U/mL Pen/Strep final (Invitrogen, Cat#: 1514022, 10000 U/mL)

5.1.2 Matrigel®, High concentration (BD Biosciences, Cat#: 354248)

IMPORTANT: All Matrigel® purchases should be submitted specifying PCR-tested LDEV-Negative Matrigel®. If not, there is a possibility of LDEV contamination which can result in LDEV+ tumors

5.2 Prepare 25% Matrigel® Working Solution

5.2.1 Chill pipettes and conical tubes in a -70°C freezer overnight and then place on wet ice prior to use.

5.2.2 Thaw Matrigel® overnight by placing a vial of Matrigel, buried in ice, in the refrigerator.

5.2.3 Using cold pipettes and tubes, make a 25% Matrigel® Working Solution with 1X F-12 Nutrient Mix + Pen/Strep, without supplementation.

5.2.4 Aliquot 25% Matrigel® either as 5-mL aliquots into chilled sterile 15-mL conical tubes or 1-mL aliquots into 2-mL sterile screw-capped tubes and place into a -20°C freezer.

6.0 PREPARATION OF MATRIGEL®-COATED FLASKS

6.1 The day before coating, remove the appropriate number of 25% Matrigel® Working Solution aliquots (prepared in Section 5.2) from the freezer and thaw overnight buried in ice, in the refrigerator.

6.1.1 The recommended volume of 2.5% Matrigel needed per well/flask are:

Plate/Flask site	Volume 2.5% Matrigel®/well or flask
96-well plate	75 µL/well
24-well plate	0.3 mL/well
6-well plate	2 mL/well
T25 flask	3.5 mL
T75 flask	6.0 mL
T162 flask	8.0 mL

6.2 Prepare 2.5% Matrigel Working solution using 25% Matrigel and 1X F-12 Nutrient Mix + Pen/Strep, without supplementation.

6.3 Coat Tissue Culture Plates/Flasks with 2.5% Matrigel® Solution

6.3.1 Matrigel® coated plates/flasks should be prepared at least 1 hour before use.

6.3.2 Using sterile procedures, coat the growth surface of the plate/flask using the recommended volume based on the size of the well/flask (Section 6.1.1). Be sure to rock the plate back and forth to completely coat the surface of the plate/flask.

6.3.3 Incubate plates/flasks for a minimum of 30 minutes at ambient temperature for polymerization. This process can be enhanced by incubating in a 37°C incubator.

6.3.4 Immediately before use, remove excess media from the flask/plate and discard taking care to not dislodge the Matrigel® coating.

- Flasks can be prepared up to 4-5 days before use and stored at 4°C following the polymerization step. Bring to ambient or 37°C before adding cells.