

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 1 of 12

TABLE OF CONTENTS

1.0	PURPOSE/SCOPE	2
2.0	SAFETY	2
3.0	RELATED DOCUMENTS	2
4.0	MATERIALS & EQUIPMENT – REVIEW ALL EQUIPMENT BASED ON METHOD OF INTEREST	2
4.1	Implant Materials & Equipment	2
4.2	Alternate Implant Methodology: Additional Materials	3
4.3	Cryopreservation Materials & Equipment	3
5.0	TISSUE PREPARATION	4
5.1	Fresh/Overnight Shipped Patient Tissue	4
5.2	Cryopreserved PDX Fragments	6
6.0	IMPLANTATION	6
6.1	Subcutaneous Implantation Procedure	6
6.2	Tumor Growth Supplements	7
7.0	COMMON PDMR ALTERNATE IMPLANT METHODOLOGIES	7
7.1	“Matrigel® Plug” Implant [1, 2]	7
7.2	Mammary Fat Pad (MFP) Implant	8
8.0	CLEAN-UP	9
9.0	EXPERIMENTAL ENDPOINTS	9
9.1	Observation	9
9.2	Tumor Assessment and Humane Euthanasia	9
9.3	Static and No-Growth Monitoring	10
10.0	CRYOPRESERVATION OF PDX MATERIAL	10
10.1	Tumor preparation for cryopreservation	10
10.2	Stepped -Rate Cryopreservation Procedure	10
10.3	Slow-rate freezing (isopropanol-based using a cryo 1°C cell-freezing container such as Mr. Frosty)	11
10.4	Slow-rate freezing (non-isopropanol based such as CoolCell)	11
	REFERENCES	12

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 2 of 12

1.0 PURPOSE/SCOPE

This Standing Operating Procedure (SOP) describes the procedures for tissue preparation and procedures for subcutaneous implantation to generate patient-derived xenografts (PDXs) and cryopreservation of material under BSL-2 safety criteria. This SOP is used/performed by the Biological Testing Branch (BTB) at NCI-Frederick, Frederick National Laboratory for Cancer Research.

2.0 SAFETY

BTB treats all patient tissue and mice carrying patient tumors as a possible health threat as the human tissue could still retain human pathogenic agents. Mice are housed in barrier facilities using full micro-isolator technique with all work, including husbandry, provided by experienced technical staff. The primary mouse strain used is the NOD.Cg-Prkdc^{scid}Il2rg^{tm1Wjl}/SzJ (NSG) which are highly susceptible to infection due to their profound immunodeficiency. All materials coming into the barrier facilities are decontaminated by autoclaving, or chemical means including the use of chlorine dioxide gas and Virkon® dip tanks for non-autoclavable items.

3.0 RELATED DOCUMENTS

Source	Title
BTB/DCTD	Fresh Tumor Collection and Handling for Generation of Patient-Derived Models

4.0 MATERIALS & EQUIPMENT – REVIEW ALL EQUIPMENT BASED ON METHOD OF INTEREST

4.1 Implant Materials & Equipment

- 4.1.1 Pipetman, 200 µL and sterile tips
- 4.1.2 Sterile scalpel and forceps
- 4.1.3 Sterile petri dishes
- 4.1.4 Matrigel® (BD BioSciences, Bedford, MA.)
 - **IMPORTANT:** All Matrigel® purchases should be submitted specifying PCR-tested LDEV-Negative Matrigel®. If not there is a possibility of LDEV contamination which can result in LDEV+ tumors as a consequence.
- 4.1.5 11-gauge Trocar (tissue implant needle), sterile (pre-chilled)
- 4.1.6 1-cc syringe with 22-gauge needle, sterile (pre-chilled)
- 4.1.7 Bone rongeur, sterilized
- 4.1.8 Rituximab (Rituxan® 10 mg/mL; Genentech)
- 4.1.9 RPMI-1640 containing primocin (working concentration: 100 µg/mL)
- 4.1.10 NOD.Cg-Prkdc^{scid}Il2rg^{tm1Wjl}/SzJ (NSG) mice, sex-matched to human patient
- 4.1.11 Isoflurane and anesthesia machine for delivery
- 4.1.12 Marcaine solution, 0.25% (Bupivacaine HCl)
- 4.1.13 Nolvasan® (disinfectant for animal surgery)

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 3 of 12

- 4.1.14** Bleach (Clorox, 5.25% Hypochlorite) diluted 1:10, 2% Virkon®, or similar disinfectant
- 4.2** Alternate Implant Methodology: Additional Materials
- 4.2.1** Estradiol pellets (Innovative Research; 3.0-mm, 90-day release, 0.18 mg)
- 4.2.2** Testosterone pellets (Innovative Research; 3.0-mm, 90-day release, 12.5 mg)
- 4.2.3** 10X PBS, sterile
- 4.2.4** Sterile water
- 4.2.5** 10 mm cloning ring, sterile (Sigma Catalog# C2059-1EA)
- 4.2.6** 6-well plate, sterile
- 4.2.7** Advanced DMEM/F12 containing 5% fetal bovine serum (FBS) and primocin (working concentration: 100 µg/mL)
- 4.2.8** Wound closure: sterile staples or sutures
- 4.3** Cryopreservation Materials & Equipment
- 4.3.1** Pipetman, 1000 µL and sterile tips
- 4.3.2** Sterile scalpel, forceps and scissors
- 4.3.3** Sterile petri dish
- 4.3.4** RPMI-1640 media containing 20% fetal bovine serum (FBS) and 10% DMSO (Burdick&Jackson Brand); made same day as use
- 4.3.5** 2-mL Cryovials and tube holder
- 4.3.6** Tumor material freshly harvested from PDX-bearing mice RPMI-1640 with primocin, maintained at 4°C by placing in sterile petri dish on wet ice
- 4.3.7** Bucket of ice (wet not dry)
- 4.3.8** Stepped Rate Cryopreservation Freezer
- An alcohol-free cell freeze container (e.g., CoolCell) or isopropanol-based (e.g., Mr. Frosty) freeze container can also be used for slow-rate freezing if a stepped rate freezer is not available.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 4 of 12

5.0 TISSUE PREPARATION

- Be sure animals are ready for implantation **immediately following** tissue preparation.
- At all patient tissue preparation steps, whenever possible some residual tissue is kept for pathology review.
- All tissue preparation should be performed in a Biological Safety Cabinet (BSC) using sterile instruments and technique. Maintain sterility at all times.

5.1 Fresh/Overnight Shipped Patient Tissue

5.1.1 Tumor Fragments and 18-g Needle Biopsies

- 5.1.1.1 Tumor tissue should be implanted as soon as practical post-collection. The tumor should be maintained at 4°C in an appropriate storage medium (e.g., CO₂ independent media containing primocin) (see PDMR shipping and collection protocol) until implantation. The PDMR has detected viable tumor material up to 96-hr post-collection using out established shipping and handling procedures, but implantation within 24 hours is desired.
- 5.1.1.2 Transfer the tumor material from the shipping/transport media into a sterile petri dish along with a small volume of transport media adequate to keep the tissue wet.
- 5.1.1.3 Using a sterile scalpel cut the tissue into 2 mm³ fragments (~30 mg). Evaluate the material prior to cutting to ensure tumor tissue is being implanted as non-tumor tissue (e.g., cartilage, muscle) can provide supporting cells but should not be implanted separately. Grossly necrotic areas are avoided.
- From resected tissue, generally 5-10 P0 mice (first passage in mouse) can be implanted. 18-gauge biopsy material is generally sufficient for 1-2 P0 mice.
 - In addition, to avoid loss of viable tumor cells, the transport media is centrifuged to pellet any shed tumor material and this is also implanted or cultured in vitro (SOP Step 5.1.2).
- 5.1.1.4 Add Rituximab to the media containing the tumor fragments to a final concentration of 1 mg/mL. Mix well and incubate 30 minutes at room temperature.
- 5.1.1.5 Place one tumor fragment into the end of a pre-chilled 11-gauge Trocar needle and add a drop of Matrigel® (~50-100 µL).
- 5.1.1.6 Keep the Trocar on ice to prevent the Matrigel® from polymerizing until implantation. polymerizing Shed tumor cells in transport media
- 5.1.1.7 The PDMR always assumes there are shed viable cells in the shipping media.
- 5.1.1.8 Transfer the transport media to a 50-mL conical-bottomed tube and centrifuge at 200-300xg for 5 minutes.
- 5.1.1.9 Disinfect the outside of the tube. Remove all but 30-50 µL of the media without disturbing the pellet – a pellet may not be visible.
- 5.1.1.10 Recap tube and gently flick to re-suspend cells in the remaining media.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 5 of 12

5.1.1.11 Using a 1-cc syringe with 22-gauge needle (both pre-chilled), draw up 100 μ L Matrigel® (kept on ice to prevent polymerizing) followed by the re-suspended cells. Avoid air bubbles in the process.

5.1.1.12 Keep syringe on ice to prevent Matrigel® from polymerizing until implantation.

5.1.2 Malignant effusions (e.g., pleural, peritoneal, etc.)

5.1.2.1 Same-day receipt of up to 1-L volumes of malignant effusions. Sterile jar containing effusion delivered packed on ice.

5.1.2.2 Use sterile technique to transfer effusion to 50-mL conical tubes or 250 ml bottles and pellet cells at 200-300xg for 5 minutes.

5.1.2.3 Remove supernatant without disturbing the cells and resuspend in the residual volume.

- If significant red blood cells are present, use hypotonic shock to lyse them as follows: (1) Fill one pipette with 9-mL of sterile water and a second pipette with 1-mL 10X PBS (or other 10X physiological salt solution). (2) Add the 9-mL sterile water to the tube containing the cells and immediately add the 1-mL of 10X salt solution – mix well.

5.1.2.4 Pool cells from all centrifuge tubes into a single conical tube and pellet the cells by centrifugation at 200xg for 5 minutes and resuspend cells in residual volume (generally 30-50 μ L). Transfer the resulting cells to a conical tip microfuge tube.

5.1.2.5 Add Rituxan stock solution to the cell suspension to a final concentration of 1 mg/mL. Incubate for 30 minutes at room temperature.

5.1.2.6 Using a 1-cc syringe with 22-gauge needle (both pre-chilled), draw up 100 μ L Matrigel® (kept on ice to prevent polymerization) followed by the re-suspended cells. Avoid air bubbles in the process.

5.1.2.7 Keep syringe on ice to prevent Matrigel® from polymerizing until implantation.

5.1.3 Bony tumor tissue (osteosarcomas, bony metastases)

5.1.3.1 Transfer the tumor material to a sterile petri dish along with some of the transport material.

5.1.3.2 Use a sterile bone rongeur and forceps to break-up the bony components of the tumor tissue.

5.1.3.3 Follow the same mouse preparation procedures as that used for SC implantation.

- Note: Bony components can be implanted without difficulty although they may not fit into a trocar, they can be inserted using forceps.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 6 of 12

5.2 Cryopreserved PDX Fragments

Important: All set-up for implantation should be in place prior to beginning to thaw the cryopreserved tissue fragments. Once the cryopreserved tissue fragments are removed from the freezer, the following steps should proceed quickly.

- 5.2.1 Remove cryovial from liquid nitrogen, transfer to lab on dry ice, and hold until ready to implant.
- 5.2.2 When ready to thaw, swirl the tube in a 37°C water bath until the contents are almost completely melted.
- 5.2.3 Spray/wipe outside of cryovial with disinfectant, then move the vial into a BSC.
- 5.2.4 Transfer PDX fragments into a petri dish containing fresh DMSO-free RPMI-1640 medium to rinse/dilute the DMSO-containing freeze media away.
- 5.2.5 Place the tumor fragment into the end of an 11-gauge Trocar needle and add a drop of Matrigel® (~50-100 µL).
- 5.2.6 Keep Trocar on ice to prevent Matrigel® from polymerizing until implantation.

6.0 IMPLANTATION

6.1 Subcutaneous Implantation Procedure

- 6.1.1 Shave the implantation site prior to anesthesia.
- 6.1.2 Anesthetize a sex-matched NSG mouse by isoflurane inhalation until non-responsive to toe-pinch. Keep mouse warm and monitor breathing during procedure.
- 6.1.3 Disinfect implantation site with Nolvasan®. We recommend to NOT use alcohol, as this can chill the animal during the procedure.
- 6.1.4 Using scissors, make a small nick incision in the inguinal area (for standard SC implants placed in the posterior axillary area) and insert the Trocar through the skin opening.
- 6.1.5 Slide the Trocar (or 1-cc syringe with 22-gauge needle) craniodorsally to about 0.5 cm caudal to the axillary area; the Trocar should slide smoothly under the skin.
 - 6.1.5.1 Notes on Implantation Using Trocar/Matrigel® Method
 - Expel the tumor fragment into the axillary area. This may require twisting the Trocar slightly to ensure the tumor fragment does not stick to the Trocar and then gently pull out the Trocar. Forceps can be used over the skin to hold the fragment in place during trocar withdrawal.
 - 6.1.5.2 Notes on Implantation Using 1-cc Syringe/Matrigel® Method
 - Inject the cell/Matrigel® inoculum into the subcutaneous space on the lateral body wall just posterior to the axillary area. Leave the needle in place in the mouse for a few seconds once injected to allow the expelled Matrigel® to solidify; this will prevent the solution from flowing out of the needle track.
- 6.1.6 The entry site will heal well and does not require a wound clip, although one may be applied if preferred.
- 6.1.7 Apply 1-2 drops 0.25% Marcaine to the nick incision to relieve pain/discomfort.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 7 of 12

6.1.8 Keep mice warm and monitor breathing until awakened from anesthesia.

6.2 Tumor Growth Supplements

6.2.1 Tumor types that may be dependent upon exogenously administered hormones (e.g., breast, prostate, ovarian, testicular) should receive subcutaneously implanted hormone pellets to aid in tumor survival and growth.

6.2.2 Estradiol pellets (Innovative Research; 3.0-mm, 90-day release, 0.18 mg) should be implanted subcutaneous into mice receiving tissue implants originating from female patients with cancer of the Breast, Ovary, Uterus, Fallopian tube, or Cervix unless contra-indicated in patient diagnoses record.

6.2.2.1 A single pellet (or 2 if desired) used in an NSG mouse should be adequate to provide the needed support for tumor growth and survival while avoiding estrogen toxicosis. Only one pellet should be used in athymic nude mice due to their sensitivity to estrogen toxicosis.

6.2.2.2 Injectable estradiol cypionate is avoided if possible as it is more likely to result in estrogen toxicosis before the tumor is able to grow.

6.2.2.3 Implant pellet at a SC site distant to the PDX implant. The pellet should be replaced every 90 days.

6.2.3 Testosterone pellets (Innovative Research; 3.0-mm, 90-day release, 12.5 mg) should be implanted subcutaneously into mice receiving tissue implants originating from male patients with cancer of the testicles, prostate, vas deferens, and seminal vesicles unless contra-indicated in patient diagnoses record.

6.2.3.1 A single pellet in an NSG mouse should be adequate to provide the needed support for tumor growth and survival.

6.2.3.2 Implant pellet at a site distant to the PDX implant. The pellet should be replaced every 90 days.

7.0 COMMON PDMR ALTERNATE IMPLANT METHODOLOGIES

PDMR implants essentially all tumor types subcutaneously (SC) per the method above except breast tumors which are implanted in the mammary fat pad (MFP). Breast cancer tumors are maintained in the MFP through all passages.

Adult Soft Tissue Sarcomas are implanted using the “Sarcoma Plug” method below as well as SC (if sufficient material) at P0. P1 tumors are then implanted at the alternate site as well as SC to determine if SC growth is supported. For adult soft tissue sarcomas, the PDMR noted a low take-rate with direct SC implantation of material. The “Sarcoma Plug” method has increased take-rate from <5% to 35%-40% across all sarcoma histologies.

Other alternate implant methodologies are used less frequently and will be posted as a separate SOP to the PDMR website.

7.1 “Sarcoma Plug” Implant [1, 2]

7.1.1 The Day Before Implantation:

7.1.1.1 Using a sterile forceps, press a sterile 10-mm cloning ring, greased side down, in the center of a well in a sterile 6-well plate to create a seal between the plate and the ring.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 8 of 12

7.1.1.2 Quickly layer ~100 μ L of Matrigel® into the cloning ring followed by addition of a small tumor fragment. Cover the fragment with an additional~ 100 μ L of Matrigel®.

7.1.1.3 Add adequate media (Advanced DMEM/F12 with 5% FBS and primocin) to fill the cloning ring and incubate overnight in a 37°C CO₂ incubator. If the ring must be incubated longer than 24-h, then the entire well should be filled with media such that the cloning ring is covered to prevent dessication.

7.1.2 Day of Implantation

7.1.2.1 Using a sterile forceps, remove the cloning ring from the 6-well plate and sterilely remove the “Sarcoma Plug” from the cloning ring. We have found using the end of a sterile 1-cc plunger works well to expel the plug from the cloning ring.

7.1.2.2 Follow the same mouse preparation procedures as that used for SC implantation with the following changes:

- Make a blunt dissection just cranial to the rear leg for access to the #4 (#9) MFP.
- Use a small spatula to “scoop” up the “Sarcoma Plug” and place it into the subcutis of a surgically anesthetized mouse, near but not within the MFP.

7.1.2.3 Close the skin and staple (or suture if preferred) the overlying skin

7.1.2.4 Apply 1-2 drops 0.25% Marcaine to the nick incision to relieve pain/discomfort.

7.1.2.5 Keep mice warm and monitor breathing until awakened from anesthesia.

7.2 Mammary Fat Pad (MFP) Implant

7.2.1 Follow the same mouse preparation procedures as that used for SC implantation with the following changes:

7.2.1.1 Make a small incision just cranial to the rear leg for visualization of the #4 (#9) MFP.

7.2.1.2 Gently lift the MFP with forceps, incise the gland crating a small pocket. Slide the Trocar into this space and expel the prepared tumor fragment into the pocket.

7.2.1.3 Return the gland to its proper position, close and staple (or suture if preferred) the overlying skin.

7.2.1.4 Apply 1-2 drops 0.25% Marcaine to the nick incision to relieve pain/discomfort.

7.2.1.5 Keep mice warm and monitor breathing until awakened from anesthesia.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 9 of 12

8.0 CLEAN-UP

- 8.1.1 All materials coming into contact with patient tissue as well as the mice carrying patient tumor samples are treated as a potential health threat (BSL-2 precautions) since the human tissues could retain human pathogenic agents even if they do not replicate in mouse cells (e.g., HIV, HPV, etc).
- 8.1.2 Flush/soak any items (e.g., tubes, syringes, petri dishes, lab mats, etc) that were in contact with human tissue with disinfectant (e.g., 10% bleach, commercial hydrogen peroxide disinfectant, 2% Virkon®) for a minimum of 10 minutes before disposal in biohazard waste or sharps containers (follow institutional guidelines and manufacturer's recommendations).
- 8.1.3 For items that can't be rinsed (e.g., micropipettors), wipe down thoroughly with bleach-soaked gauze or other appropriate disinfectants.

9.0 EXPERIMENTAL ENDPOINTS

9.1 Observation

- 9.1.1 Mice are weighed 1-3 times per week as a gauge of animal health and more frequently if clinical signs are noted.
- 9.1.2 Mice are observed at least once daily and tumor growth is monitored with caliper measurements at least once weekly. The tumor weight (mg) is calculated using the following formula [3]:

$$Tumor\ weight\ (mg) = \frac{(tumor\ length) \times (tumor\ width)^2}{2}$$

- 9.1.3 Potential growth of P0 (first passage in mouse) tumors at NCI is monitored for up to 300 days after implantation. Mice not showing progressive tumor growth at 300 days post-implant are euthanized.
- 9.1.4 Antemortem determination for harvesting tumor tissues is based upon tumor size as well as clinical signs. Animals without palpable tumors are evaluated clinically and by monitoring body weights.

9.2 Tumor Assessment and Humane Euthanasia

- 9.2.1 Tumors are harvested immediately following euthanasia or under terminal isoflurane anesthesia.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 10 of 12

- 9.2.2** Palpable tumors are evaluated 2-3 times per week and efforts are made to identify animals that may experience adverse health effects before the next observation (and euthanize accordingly).
- 9.2.2.1 A superficial tumor burden of 4000 mg is an endpoint; however, tumor harvests in the 1000-2000 mg size range are preferred as there is typically less necrotic tissue present.
- 9.2.2.2 An animal with a tumor that grows in a manner that limits their mobility, impedes access to food and water, or has a tumor that is showing the onset of skin necrosis is euthanized.
- 9.2.2.3 Tumors may be red to purple in color when highly vascularized and this appearance does not warrant euthanasia. When the skin overlying the tumor shows evidence of necrosis or other compromise, euthanasia is warranted.
- 9.2.2.4 Animals with tumors arising at sites other than the implantation site will be euthanized based mainly upon gross observations such as weight loss or gain, lethargy, respiratory distress, inactivity, ruffled hair coat, decreased body temperature, inappetance, etc.
- 9.2.2.5 Some internal tumors (prostate, spleen, renal) may be palpable and will be checked regularly so that the tumor does not become large and debilitating.

9.3 Static and No-Growth Monitoring

- 9.3.1** Animals with no signs of tumor growth are euthanized at 300 days post-implant.
- 9.3.2** NCI has observed numerous tumors that we have categorized as “Static” tumors. These tumors reach a steady-state tumor burden and appear to stop growing in size; possibly due to walling off of the tumor in the animal. These mice are euthanized and fragments of the tumor are passaged to determine if tumor growth will re-initiate.

10.0 CRYOPRESERVATION OF PDX MATERIAL

10.1 Tumor preparation for cryopreservation.

- 10.1.1** In a BSC, harvest the tumor into a sterile petri dish containing a small volume of RPMI-1640 with primocin.
- 10.1.2** Using a sterile scalpel, cut the tissue into ~2-3 mm³ fragments (~30mg).
- 10.1.3** With the cryovials in a tube rack held on ice, add 3-5 fragments to each tube then add 1-mL of freeze medium (RPMI 1640/5% FBS/10% DMSO) to each tube without exceeding the tube fill volume.
- 10.1.4** Apply the cap to the cryovials, seal well. Wipe the exterior with disinfectant then place into wet ice until ready to begin stepped--rate freezing.

10.2 Stepped -Rate Cryopreservation Procedure

- 10.2.1** Mechanically controlled stepped-rate freeze cryopreservation is recommended, when available, as it is believed to result in a lower loss of cell viability due to the decreased formation of ice crystals.
- 10.2.2** Always follow the manufacture’s guidelines for operation. General stepped-rate cryopreservation parameters used at by the PDMR are:

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 11 of 12

- Decrease 1°C/minute down to -4°C
- Decrease 25°C/minute down to -40°C
- Increase 15°C/minute up to -12°C
- Decrease 1°C/minute down to -40°C
- Decrease 10°C/minute down to -90°C

10.2.3 Transfer the cryopreserved vials to liquid nitrogen storage.

10.3 Slow-rate freezing (isopropanol-based using a cryo 1°C cell-freezing container such as Mr. Frosty)

10.3.1 Follow the manufacturer's instructions as provided for the specific cryopreservation device.

- PDX material should be held on wet ice in cryovial tubes until ready for placement into the slow-rate freeze container.
- The base of the cryo-container is filled with isopropanol per the manufacturer's recommendation and the tube holder is placed on top.
- Transfer the cryovials filled with fragments/freeze media from the ice-bucket into the tube holder of the cryo-container, screw the lid securely onto the cryo-container, and place at -80°C for a minimum of 4 hrs, but most commonly for overnight.
 - If the tubes are held on ice, the cryo-container can be at refrigerator temperatures. If the cryovials are processed at room temperature rather than 4°C, then the cryo-container should be at room temperature.
- Vials should be transferred to the vapor phase of a liquid nitrogen tank as soon as practical after the 4-hr freeze step, typically the following morning. In no case, should the vials be held longer than 3 days at -80°C before transfer into the vapor phase of a liquid nitrogen storage tank.

10.4 Slow-rate freezing (non-isopropanol based such as CoolCell)

10.4.1 Follow the manufacturer's instructions as provided for the specific cryopreservation device.

Implantation and Cryopreservation of Tissue for PDX Generation		
Laboratory:	Biological Testing Branch	
Effective Date:	4/17/2017	Page 12 of 12

REFERENCES

1. Kenney, S., et al., *ASPS-1, a novel cell line manifesting key features of alveolar soft part sarcoma*. J Pediatr Hematol Oncol, 2011. **33**(5): p. 360-8.
2. Vistica, D.T., et al., *Therapeutic vulnerability of an in vivo model of alveolar soft part sarcoma (ASPS) to antiangiogenic therapy*. J Pediatr Hematol Oncol, 2009. **31**(8): p. 561-70.
3. Plowman J, D.D., Hollingshead M, et al., *Human tumor models in NCI drug development*. Anticancer Drug Development Guide: Preclinical Screening, Clinical Trials, and Approval. 1997.